

Eldex

Optos Series HOT HEAD

Quick Heating, Easy Removal Heating Accessory for Eldex Optos Pumps



Shown with optional controller

1. Easy On/Off design for access to head and valves for maintenance.
2. 120V and 240V AC models for simple implementation worldwide.
3. Integrated bi-metal thermostat for safety shutoff.
4. Able to install using our optional Heated Head PID Controller or connect directly to your own PLC controls.

- The HOT HEAD is made out of a PEEK body with an aluminum insert for quick heating and durability.
- Quick release latch for easy removal of heater from head for easy maintenance.
- PID controller available for easily setting and maintaining a constant temperature on the head assembly.

Eldex Corporation
30 Executive Ct.
Napa, CA 94558-6278
Tel: (707) 224.8800 | (800) 969.3533
Fax: (707) 224.0688
Email: sales@eldex.com

Temperature	Up to 250 F (121C)
Hot Head Voltage	120VAC or 240VAC
Current	0.8A
Material	PEEK with Aluminum Insert
Safety Switch	Off @ 300F; ON @ 250F
Thermocouple	K-Style
Controller Voltage Input	120–240 VAC